




..... **Laser GWYYb**

Medical Laser is accessed outside of the exam except when using the Retina Exam.
In the Retina Exam, you can access the Medical Laser screen with **LP** on the right side of screen.



This opens the same Medical Laser Screen that is accessed from the Clock Screen.
The Medical Laser Screen requires all of the entries to be made on every procedure.
None of the previous information from past procedures can be carried forward to the new exam.



PROCEDURE NOTES Copyright 1999-2010 Medflow, Inc.

Patient Details: Diane Wilton, 20, 12/20/2011

DIAGNOSES: 370.22 MACULAR KERATITIS

Medications Administered: Mydracyl BI (No 25%) Pilo 2%

Laser Type: ARGON

Diode Used: 809nm (809nm continuous), 810nm (micropulse)

Application: Spots

Number of Applications: 1, Duration: 1Min, Spot Size: 1mm

PostOpMeds: Acular opth sol, 1 gr Q4h v/a to operative eye

Discharge Status: To Home, Alert & Oriented X3

Follow Up: Follow up in 24 hours, Return in 4 days, Return in 1 month, Return in 3 weeks, Return in 3 months, Return in A.M.

Buttons: View Patient Documents, View Images, Doctor Screen, Schedule, CANCEL, DRIV Fluidat, DRIV Gonio, DRIV Cones, Workup, PN, FDC, EVALUATION TIME (min), SAVE/EXIT

Choose the laser procedure from the list and then select **OK**. This will open the Laser Screen.
 From this screen it is possible to load past information into today's procedure and then update the information for today's visit.

If previous lasers had been performed on this patient, you can select EDIT EXISTING LASER EXAMS

Double click on that entry, fill in the information on today's screen and make whatever changes necessary. Pick the laser type in the lower left hand corner
 You can make any notes in the IntraOpNote area by free typing. Then you would **Save and Generate Document**.



This is what the document will look like once generated. The provider can **SAVE/SIGN/EXIT**, the technician will **SAVE/EXIT** and this will send the document to the provider's sign chart box.

Medflow Document Editor Copyright: 1999-2009 Medflow, Inc.

TEMPLATES Refresh Insert Header Info Wilson, Diane GENDER:F DOB:4/30/1965

a-laser-Triab.rtf Argon-laser-Treb.rtf BSCAN.RTF CATARACT.rtf CATARACTOPNOTE.RTF

SAVE/SIGN/EXIT SAVE DOCUMENT RECIPIENT LIST 0 recipients

set printer Black Courier New 10

Medflow, Inc
6739 A Fairview Rd
Charlotte, NC 28212
704-927-9800 704-927-1234
12/20/2011

PATIENT NAME
Wilson, Diane GENDER:F DOB:4/30/1965

| ADDRESS | CITY | STATE | ZIP |
|---------------------|-----------|-------|-------|
| 6739A Fairview Road | Charlotte | NC | 28210 |

ACCOUNT#: 20
DIAGNOSIS: 366.16 CATARACT
NUCLEAR SCLEROSIS OU
PROCEDURE: 66821 YAG CAPSULOTOMY
HEIGHT: 0 inches
WEIGHT: 0 lbs
EYE: OD
TREATMENT SPOT SIZE: 12 microns
TREATMENT POWER: 200 mW/cm²
TREATMENT DURATION: 20 seconds
NUMBER OF SPOTS: 52

Doctor: Kevin Cornea M.D.

Prepared by: EMP

PRINT FAX CANCEL SPELL CHECK SELECTED TEXT SAVE/EXIT

envelope